

All NEW Diamond® Midi™

5mL Centrifuge Tubes

An efficient and economical alternative to 15mL conical tubes when working with sample volumes between 2.0mL and 5.0mL.

With the same diameter and conical shape as 15mL tubes, our shorter 5.0mL Diamond® Midi™ centrifuge tubes:

- Reduce total cost to the lab for purchase and use
- Reduce environmental impact during production and disposal
- Require no special equipment or racks
- Reduce pipet shaft contamination since there is no need to insert the pipettor barrel into the tube to retrieve samples
- Available non-sterile or sterile and are RNase, DNase and pyrogen free
- Spin up to 25,000 x g when used with proper support in appropriate rotor

Diamond® Midi™ 5mL Screw Cap Centrifuge Tubes

- Single turn threaded caps seal easily and securely
- Graduated in 0.5mL and 0.2mL increments
- Available sterile in bags or recyclable paperboard racks and non-sterile with tubes and caps packed separately
- Suitable for use from -80° to 121°C for low temperature sample storage, general laboratory uses, or boiling procedures
- Polypropylene (PP) with polypropylene (PP) screw cap

Diamond® Midi™ 5mL Snap Cap Centrifuge Tubes

- Provides the convenience of a snap cap microcentrifuge tube in a tube that holds larger volumes
- Attached snap cap facilitates one handed operation
- Pierceable cap allows easy access by needle or syringe
- Dual graduation scales in increments of 0.5mL and 0.2mL enable more precise measurement when required
- Suitable for use from -80° to 121°C for low temperature sample storage and general laboratory use
- Natural tubes include boil proof closure
- Seven color options
- Polypropylene (PP)

**Buy 3 Get 1
FREE OFFER***



111578



111580



111580B

***HURRY! Limited time introductory offer!**

Order through your authorized Globe Scientific dealer between October 9th and December 31st, 2019. Limit 4 free units per lab. Mix & match. Free items of equal or lesser value.

See redemption form on next page.



www.globescientific.com

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

Diamond® Midi™ 5mL Microcentrifuge Tube REDEMPTION FORM



Buy any 3 units of Diamond® Midi™ Centrifuge Tubes and get a fourth unit **FREE!***

Please follow the instructions below to redeem promotions.

1. Order 3 units of Diamond® Midi™ centrifuge tubes between October 9 and December 31, 2019.
2. Complete the form below, then fax this redemption form and proof-of-purchase (dealer invoice or packing slip) to Globe Scientific, 201-599-1406 or e-mail to promotions@globescientific.com

Please allow up to 4 weeks for delivery of free products. All redemption forms must be received no later than 30 days after the promotion's expiration date. Offer valid for shipment within the contiguous U.S. only.

Enter Account Information

Name _____ Date _____

Organization _____ Dealer _____

Ship to Address _____

Room/Bldg./Dept. _____

City _____ State _____ Zip _____

Phone (required) _____ Email (required) _____

Enter Offer Information

Item No. of Purchased Item(s)	Qty. Purchased	Item No. of Requested Eligible FREE Item(s)	Quantity FREE

*Mix & Match. Free item of equal or lesser value. Limit 4 free units per lab. Attach additional sheet if necessary. Promo code **Midi**.

Proof of Purchase (authorized dealer invoice or packing slip) MUST accompany this redemption form!

Redemption of this offer allows Globe Scientific to periodically e-mail user about new products and promotions from Globe Scientific. Globe Scientific will never share or distribute your contact information with third parties.